

Product Datasheet

PPA1 Antibody (orb1517342)

| | |
|----------------------------|---------------------------------------------------------------|
| Description | PPA1 Monoclonal Antibody |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FC, IF, WB |
| Immunogen | Recombinant Human Inorganic pyrophosphatase protein (1-289AA) |
| Target | PPA1 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Uniprot ID | Q15181 |
| Expiration Date | 12 months from date of receipt. |