

Product Datasheet

Myc-tag Antibody (orb1517357)

Description	Myc-tag Antibody
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Myc synthetic peptide conjugated to KLH.
Target	Myc-tag
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.05% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.