

Product Datasheet

HA-tag Antibody (orb1517359)

Description	HA-tag Antibody
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	HA synthetic peptide conjugated to KLH.
Target	HA-tag
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.05% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.