

Product Datasheet

FLAG-tag Antibody (orb1517360)

| | |
|----------------------------|--|
| Description | FLAG-tag Antibody |
| Species/Host | Rabbit |
| Reactivity | All |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | DDDDK synthetic peptide conjugated to KLH. |
| Target | FLAG-tag |
| Preservatives | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.05% sodium azide. |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Expiration Date | 12 months from date of receipt. |