

Product Datasheet

Rubisco Antibody (orb1517389)

Description	Rubisco Antibody
Species/Host	Mouse
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Target	Rubisco
Preservatives	pH 7.4 150mM NaCl 0.02% sodium azide and 50% glycerol.a
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.