

Product Datasheet

FLAG tag Antibody (orb1517391)

Description	FLAG tag Antibody
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide: DYKDDDDK conjugated to KLH
Target	FLAG tag
Preservatives	Rabbit IgG in pH7.3 PBS, 0.05% NaN ₃ , 50% Glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.