

Product Datasheet

AmCyan Antibody (orb1517395)

Description	AmCyan Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	The AmCyan Mouse Polyclonal Antibody detects AmCyan & AmCyan Tag protein.
Target	AmCyan
Preservatives	pH 7.4 150mM NaCl 0.02% sodium azide and 50% glycerol.al
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.