

Product Datasheet

MCP Antibody (orb1517402)

Description	MCP Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human cytomegalovirus (strain AD169) (HHV-5) (HCMV) Major capsid protein (862-1048aa)
Target	MCP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16729
Expiration Date	12 months from date of receipt.