

Product Datasheet

Non-structural polyprotein Antibody (orb1517406)

| | |
|----------------------------|---|
| Description | Non-structural polyprotein Antibody |
| Species/Host | Rabbit |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant Venezuelan equine encephalitis virus (strain Trinidad donkey) (VEEV) Non-structural polyprotein (1-288aa) |
| Target | Non-structural polyprotein |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P27282 |
| Expiration Date | 12 months from date of receipt. |