

Product Datasheet

Polyprotein Antibody (orb1517414)

Description	Polyprotein Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Ugandan cassava brown streak virus Polyprotein (1407-1423aa)
Target	Polyprotein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	C7C699
Expiration Date	12 months from date of receipt.