

Product Datasheet

LBA Antibody (orb1517416)

Description	LBA Antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Glycine max (Soybean) (Glycine hispida) Leghemoglobin A (10-31aa)
Target	LBA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02238
Expiration Date	12 months from date of receipt.