

Product Datasheet

mpt64 Antibody (HRP) (orb1517419)

| | |
|----------------------------|---|
| Description | mpt64 Antibody, HRP conjugated |
| Species/Host | Rabbit |
| Reactivity | Bacteria |
| Conjugation | HRP |
| Tested Applications | ELISA |
| Immunogen | Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64 protein (24-228AA) |
| Target | mpt64 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P9WIN9 |
| Expiration Date | 12 months from date of receipt. |