

Product Datasheet

ORF3 Antibody (orb1517425)

| | |
|----------------------------|---|
| Description | ORF3 Antibody |
| Species/Host | Rabbit |
| Reactivity | Drosophila |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant Drosophila melanogaster (Fruit fly) RNA-directed DNA polymerase (666-895aa) |
| Target | ORF3 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q5BN19 |
| Expiration Date | 12 months from date of receipt. |