

Product Datasheet

A29L Antibody (orb1517430)

Description	A29L Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Monkeypox virus (strain Zaire-96-I-16) (MPX) A29L protein
Target	A29L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q77HM6
Expiration Date	12 months from date of receipt.