

Product Datasheet

GUCY2C Antibody (orb1517458)

Description	GUCY2C Monoclonal Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Recombinant Human GUCY2C protein
Target	GUCY2C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	mIgG2a
Clonality	Monoclonal
Uniprot ID	P25092
Dilution Range	Recommended dilution: FC:1:50-1:200
Expiration Date	12 months from date of receipt.