

Product Datasheet

TNFRSF9/4-1BB Antibody (orb1517471)

Description	TNFRSF9/4-1BB Monoclonal Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human TNFRSF9 protein
Target	TNFRSF9/4-1BB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	hIgG2
Clonality	Monoclonal
Uniprot ID	Q07011
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.