

## Product Datasheet

### HA Tag Antibody (orb1517495)

<b>Description</b>	HA Tag Antibody
<b>Species/Host</b>	Goat
<b>Conjugation</b>	Agarose
<b>Tested Applications</b>	IP
<b>Target</b>	HA Tag
<b>Preservatives</b>	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
<b>Form/Appearance</b>	Whole IgG
<b>Storage</b>	2 - 8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Format: Whole IgG
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	IP - 15 - 25 µl of gel slurry per 0.1 to 1 mg of protein lysate or extract
<b>Expiration Date</b>	12 months from date of receipt.