

Product Datasheet

ECS (DYKDDDDK) Tag Antibody (orb1517499)

Description	ECS (DYKDDDDK) Tag Antibody
Species/Host	Goat
Conjugation	Agarose
Tested Applications	IP
Target	ECS (DYKDDDDK) Tag
Preservatives	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	IP - 15 - 25 µl of gel slurry per 0.1 to 1 mg of protein lysate or extract
Expiration Date	12 months from date of receipt.