

Product Datasheet

PRKAR2A/PKA-RIIalpha Antibody (orb1517804)

Description	PRKAR2A/PKA-RIIalpha Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Between 50 and 100
Target	PRKAR2A/PKA-RIIalpha
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P13861
NCBI	NP_004148.1
Dilution Range	IHC - 1:100 - 1:500
Expiration Date	12 months from date of receipt.