

Product Datasheet

ERCC5/XPG Antibody (orb1517922)

| | |
|----------------------------|--|
| Description | ERCC5/XPG Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC |
| Immunogen | Between 650 and 700 |
| Target | ERCC5/XPG |
| Preservatives | Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide |
| Form/Appearance | Whole IgG |
| Concentration | 200 µg/ml |
| Storage | 2 - 8°C |
| Note | For research use only |
| Application notes | Format: Whole IgG |
| Clonality | Polyclonal |
| Uniprot ID | P28715 |
| NCBI | NP_000114.2 |
| Dilution Range | IHC - 1:100 - 1:500 |
| Expiration Date | 12 months from date of receipt. |