

Product Datasheet

METTL3/MT-A70 Antibody (orb1518225)

Description	METTL3/MT-A70 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Between 1 and 50
Target	METTL3/MT-A70
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q86U44
NCBI	NP_062826.2
Dilution Range	IHC - 1:100 - 1:500
Expiration Date	12 months from date of receipt.