

Product Datasheet

EFTUD2/SNRP116 Antibody (orb1518550)

Description	EFTUD2/SNRP116 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Between 1 and 50
Target	EFTUD2/SNRP116
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q15029
NCBI	NP_004238.2
Dilution Range	IHC - 1:100 - 1:500
Expiration Date	12 months from date of receipt.