

Product Datasheet

SMAR1/BANP Antibody (orb1518589)

Description	SMAR1/BANP Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Between 469 and C-terminus
Target	SMAR1/BANP
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	50 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8N9N5
NCBI	Q96GJ7
Dilution Range	IHC - 1:100 - 1:500
Expiration Date	12 months from date of receipt.