

Product Datasheet

XRCC1, Phospho (S518/T519/T523) Antibody (orb1518638)

Description	XRCC1, Phospho (S518/T519/T523) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Surrounding serine 518 and threonine 519/523
Target	XRCC1, Phospho (S518/T519/T523)
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	50 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P18887
NCBI	NP_006288.2
Dilution Range	IHC - 1:100 - 1:500; IHC-IF - 1:50 - 1:500
Expiration Date	12 months from date of receipt.