

Product Datasheet

Aurora A, Phospho (T288) Antibody (orb1518717)

Description	Aurora A, Phospho (T288) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC
Immunogen	Surrounding threonine 288
Target	Aurora A, Phospho (T288)
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	250 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O14965
NCBI	NP_003591.2
Dilution Range	ICC-IF - 1:100 - 1:500
Expiration Date	12 months from date of receipt.