

Product Datasheet

SMC1, Phospho (S966) Antibody (orb1518810)

Description	SMC1, Phospho (S966) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Surrounding serine 966
Target	SMC1, Phospho (S966)
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	250 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q14683
NCBI	NP_006297.2
Dilution Range	IHC - 1:100 - 1:500
Expiration Date	12 months from date of receipt.