

Product Datasheet

UBQLN4/CIP75/Ubiquilin 4 Antibody (BSA Free) (orb1519631)

Description	UBQLN4/CIP75/Ubiquilin 4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	Residues 2-601 (full length)
Target	UBQLN4/CIP75/Ubiquilin 4
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, rAlbumin-Free
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	Q9NRR5
NCBI	NP_064516.2
Dilution Range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , CHIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF:
Expiration Date	12 months from date of receipt.