


Product Datasheet

CDK1, Phospho (Y15) Antibody (orb1519743)

Description	CDK1, Phospho (Y15) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Surrounding tyrosine 15
Target	CDK1, Phospho (Y15)
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	40 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P06493
NCBI	NP_001777.1
Dilution Range	WB - 1:1,000; IP - 20 µl/1 mg lysate
Expiration Date	12 months from date of receipt.



Detection of human Phospho CDK1 (Y15) by...