

Product Datasheet

CEACAM1/5 Antibody (BSA Free) (orb1519920)

| | |
|----------------------------|--|
| Description | CEACAM1/5 Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IHC, IP, WB |
| Immunogen | Residues 35-433 (ECD) |
| Target | CEACAM1/5 |
| Preservatives | Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, rAlbumin-Free |
| Form/Appearance | Whole IgG |
| Concentration | 1000 µg/ml |
| Storage | 2 - 8°C |
| Note | For research use only |
| Application notes | Format: Whole IgG |
| Clonality | Recombinant |
| Uniprot ID | P13688 , P06731 |
| NCBI | NP_001020083.1 |
| Dilution Range | WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF: |
| Expiration Date | 12 months from date of receipt. |