



Product Datasheet

53BP1 Antibody (orb1519978)


Description	53BP1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Immunogen	Between 350 and 400
Target	53BP1
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	Q12888
NCBI	NP_005648.1
Dilution Range	WB - 1:1000; IP - 20 µl/mg lysate; IHC - 1:100 to 1:500. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human 53BP1 by western blot...



Detection of human 53BP1 by immunohistoc...



Detection of mouse 53BP1 by immunohistoc...