

Product Datasheet

AHA1 antibody (orb152015)

Description	Rabbit polyclonal to Aha1 (Biotin). Aha1 is a member of the Hsp90 cochaperone family, and is thought...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC, IP, WB
Immunogen	Recombinant Full Length Mouse Aha1 Protein
Target	AHA1
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µl/ml of SPC-183 was sufficient for detection of Aha1 in 10 µg of mixed human cell lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	38kDa
Uniprot ID	Q8BK64
NCBI	NP_666148.1
Entrez	217737
Dilution Range	WB (1:1000), ICC/IF (1:60)
Expiration Date	12 months from date of receipt.