



Product Datasheet

GMIP Antibody (orb1520228)


Description	GMIP Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 920 and C-terminus
Target	GMIP
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9P107
NCBI	NP_057657.2
Dilution Range	WB - 1:10,000 - 1:25,000; IP - 20-40 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human GMIP by western blot....



Detection of mouse GMIP by western blot....



Detection of human GMIP by western blot ...