

Product Datasheet

HSP60 (*P. falciparum*) antibody (orb152048)

Description	Rabbit polyclonal to Hsp60 (APC). In both prokaryotic and eukaryotic cells, the misfolding and aggre...
Species/Host	Rabbit
Reactivity	Bacteria, <i>E. coli</i>
Conjugation	APC
Tested Applications	ICC, IHC, WB
Immunogen	Recombinant full length PfHSP60
Target	HSP60
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1.83 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	0.9 µg/ml of SPC-185 was sufficient for detection of PfHSP60 in 20 µg of <i>P. falciparum</i> lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	60kDa
Uniprot ID	P34940
NCBI	XM_001347402.1
Dilution Range	WB (1:2000)
Expiration Date	12 months from date of receipt.