

## Product Datasheet

### HSP60 (*P. falciparum*) antibody (orb152054)

<b>Description</b>	Rabbit polyclonal to Hsp60 (RPE). In both prokaryotic and eukaryotic cells, the misfolding and aggre...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Bacteria, <i>E. coli</i>
<b>Conjugation</b>	RPE
<b>Tested Applications</b>	ICC, IHC, WB
<b>Immunogen</b>	Recombinant full length PfHSP60
<b>Target</b>	HSP60
<b>Preservatives</b>	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
<b>Concentration</b>	1.83 mg/ml
<b>Storage</b>	Conjugated antibodies should be stored at 4°C
<b>Note</b>	For research use only
<b>Application notes</b>	0.9 µg/ml of SPC-185 was sufficient for detection of PfHSP60 in 20 µg of <i>P. falciparum</i> lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
<b>Clonality</b>	Polyclonal
<b>MW</b>	60kDa
<b>Uniprot ID</b>	<a href="#">P34940</a>
<b>NCBI</b>	<a href="#">XM_001347402.1</a>
<b>Dilution Range</b>	WB (1:2000)
<b>Expiration Date</b>	12 months from date of receipt.