


Product Datasheet

MCTP2 Antibody (orb1520604)

Description	MCTP2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 50 and 100
Target	MCTP2
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q6DN12
NCBI	Q6DN12.3
Dilution Range	IP - 8 - 15 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human MCTP2 by western blot...