


## Product Datasheet

### GORAB Antibody (orb1520609)

<b>Description</b>	GORAB Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IP
<b>Immunogen</b>	Between 319 and C-terminus
<b>Target</b>	GORAB
<b>Preservatives</b>	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
<b>Form/Appearance</b>	Whole IgG
<b>Concentration</b>	1000 µg/ml
<b>Storage</b>	2 - 8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Format: Whole IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q5T7V8</a>
<b>NCBI</b>	<a href="#">Q5T7V8.1</a>
<b>Dilution Range</b>	IP - 8 - 15 µl/mg lysate
<b>Expiration Date</b>	12 months from date of receipt.



Detection of human GORAB by western blot...