



## Product Datasheet

### GLOD4 Antibody (orb1520663)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | GLOD4 Antibody   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse   |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | IP, WB   |
| <b>Immunogen</b>           | Between 263 and C-terminus   |
| <b>Target</b>              | GLOD4  |
| <b>Preservatives</b>       | Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide |
| <b>Form/Appearance</b>     | Whole IgG  |
| <b>Concentration</b>       | 1000 µg/ml   |
| <b>Storage</b>             | 2 - 8°C  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Format: Whole IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <a href="#">Q9HC38</a>   |
| <b>NCBI</b>                | <a href="#">Q9HC38.1</a>   |
| <b>Dilution Range</b>      | WB - 1:2,000 - 1:10,000; IP - 20-40 µl/mg lysate                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |



Detection of human and mouse GLOD4 by western blot...



Detection of human GLOD4 by western blot...