



Product Datasheet

UAP1 Antibody (orb1520700)

Description	UAP1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 1 and 50
Target	UAP1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q16222
NCBI	Q16222.3
Dilution Range	WB - 1:1,000 - 1:5,000; IP - 3 - 6 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human UAP1 by western blot....



Detection of human UAP1 by western blot ...