



Product Datasheet

CLGN Antibody (orb1520719)

Description	CLGN Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 1 and 50
Target	CLGN
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O14967
NCBI	O14967.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 8 - 15 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human CLGN by western blot....



Detection of human CLGN by western blot ...