



Product Datasheet

IP6K1 Antibody (orb1520726)

Description	IP6K1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 75 and 125
Target	IP6K1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q92551
NCBI	Q92551.3
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 8 - 15 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human IP6K1 by western blot...



Detection of human IP6K1 by western blot...