



Product Datasheet

NCLN Antibody (orb1520736)


Description	NCLN Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 375 and 425
Target	NCLN
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q969V3
NCBI	Q969V3.2
Dilution Range	WB - 1:1000; IP - The antibody contained within orb1520736 has been qualified for use in immunoprecipitation.
Expiration Date	12 months from date of receipt.



Detection of human NCLN by western blot....



Detection of mouse NCLN by western blot....



Detection of human NCLN by western blot ...