


Product Datasheet

FAHD1 Antibody (orb1520746)

Description	FAHD1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 171 and C-terminus
Target	FAHD1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q6P587
NCBI	Q6P587.2
Dilution Range	IP - 8 - 15 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human FAHD1 by western blot...