

Product Datasheet

ATP50 Antibody (orb1521116)

Description	ATP50 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 1 and 50
Target	ATP50
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P48047
NCBI	P48047.1
Dilution Range	WB - Not recommended. Use rabbit anti-ATP50 antibody A305-419A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.