

Product Datasheet

VBP1 Antibody (orb1521146)

Description	VBP1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 147 and C-terminus
Target	VBP1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P61758
NCBI	P61758.3
Dilution Range	WB - Not recommended. A305-403A has not performed satisfactorily when used for WB of VBP1 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of VBP1.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.