

Product Datasheet

HSP110 antibody (orb152134)

Description

Rabbit polyclonal to Hsp110 (Biotin). Hsp110

known as t...

Species/Host

Rabbit

Reactivity

Bovine, Fish, Hamster, Human, Monkey, Mouse, Rat, Sheep, Yeast

Conjugation

Biotin

Tested Applications

ICC, IHC, WB

Immunogen

Synthetic peptide derived from the sequence of hamster HSP110; sequence identical to human and mouse

Target

HSP110

Preservatives

 PBS pH 7.2, 50% glycerol, 0.09% sodium azide
 *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SPC-195 was sufficient for detection of HSP110 in 10 µg of human cell line mixed lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

110kDa

Uniprot ID
[Q92599](#)
NCBI
[NP_006635.3](#)
Entrez
[10808](#)
Dilution Range

WB (1:1000), IHC (1:100)

Expiration Date

12 months from date of receipt.