

Product Datasheet

Codanin 1/CDAN1 Antibody (orb1522010)

Description	Codanin 1/CDAN1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 1075 and 1125
Target	Codanin 1/CDAN1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8IWY9
NCBI	Q8IWY9.4
Dilution Range	WB - Not recommended. Use rabbit anti-Codanin 1/CDAN1 antibody A304-952A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.