
Product Datasheet

SH3KBP1 Binding Protein 1/SHKBP1 Antibody (orb1522127)

Description

SH3KBP1 Binding Protein 1/SHKBP1

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 75 and 125
Target	SH3KBP1 Binding Protein 1/SHKBP1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8TBC3
NCBI	Q8TBC3.2
Dilution Range	WB - Not recommended. orb1522127 has not performed satisfactorily when used for WB of SH3KBP1 Binding Protein 1/SHKBP1 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of SH3KBP1 Bi.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.