



Product Datasheet

R3HDM1/R3HDM Antibody (orb1522373)

Description	R3HDM1/R3HDM Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 500 and 550
Target	R3HDM1/R3HDM
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q15032
NCBI	NP_056176.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human and mouse R3HDM1/R3HDM...



Detection of human R3HDM1/R3HDM by western blot...