



Product Datasheet

CCT7 Antibody (orb1522445)

Description	CCT7 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	Between 1 and 50
Target	CCT7
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q99832
NCBI	NP_006420.1
Dilution Range	WB - Not recommended. Use rabbit anti-CCT7 antibody A304-730A.; IP - 2 - 10 µg/mg lysate; IHC - 1:2,000 - 1:10,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human CCT7 by western blot ...



Detection of human CCT7 by immunohistochem...