


Product Datasheet

C9orf78 Antibody (orb1522668)

Description	C9orf78 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 239 and C-terminus
Target	C9orf78
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9NZ63
NCBI	NP_057604.1
Dilution Range	WB - Not recommended. A304-613A has not performed satisfactorily when used for WB of C9orf78 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of C9orf78.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human C9orf78 by western bl...